

Denaturation	95 °C	30 s	35
Annealing/extension/data acquisition	60 °C	45 sec	35
Melt curve analysis	60-95 °C		

34l Add the appropriate DNA copies for each standard when setting up the qPCR plate in the qPCR program, as follows:

Standard	dsDNA molecules/μl
Standard 1	1.2×10^7
Standard 2	1.2×10^6
Standard 3	1.2×10^5
Standard 4	1.2×10^4
Standard 5	1.2×10^3
Standard 6	1.2×10^2

35l Analyze qPCR results. Multiply the average of duplicate values by the dilution factor and by the five-fold dilution factor of the qPCR reaction, as follows: Total copies/ μ l = # * dilution factor.

36l Pool library for MiSeq. Pool all the samples in one library at equimolar concentrations. 1X pooled library should be in a total volume of 5 μ l, $\sim 8 \times 10^9$ molecules.

37l Denature the pooled library ($\sim 8 \times 10^9$ molecules) by adding 5 μ l of NaOH 0.2N and incubate at room temperature for 5 min. Then, add 940 μ l of Hyb buffer (supplied with MiSeq® Reagent Kit v3 (600 cycle)).

38l Prepare the Phix control V3 (PhiX Control V3 KIT) as follows: mix 2 μ l of 10 nM PhiX control with 3 μ l of Tris-HCl 10 mM + 0.1% Tween-20, denature with 5 μ l of NaOH 0.2N and incubate at room temperature for 5 min. Add 990 μ l of Hyb buffer, to generate 20 pM PhiX. Then, make a 12.5 pM PhiX dilution, by mixing 375 μ l of the 20 pM PhiX with 225 μ l of Hyb buffer. Add 100 μ l of the 12.5 pM Phix to the denatured library.

39l Clean the Flow Cell (supplied with MiSeq® Reagent Kit v3 (600 cycle)) with ultra-pure water, dry with lens tissues, followed by cleaning with alcohol wipes and lens tissue.

40l Load and sequence library using a MiSeq 600-cycle v3 kit according to manufacturer's instructions using MiSeq system. Sequencing is performed with 150 bp paired-end reads and 8 bp dual-index reads.

41l After sequencing, copy the demultiplexed output FASTQ files to a location accessible to CIRCLE-seq/CHANGE-seq analysis pipeline.